



## CircuFlow™ 5150 Intermittent Pneumatic Compression (IPC) Pump

The CircuFlow™ 5150 Intermittent Pneumatic Compression (IPC) Pump is used to treat chronic lymphatic and venous disorders of the lower extremities and chronic lymphatic disorders of the upper extremities.



### Key Features:

- FDA and CE Approved
- Use indicator timer which shows patient use and compliance
- Quiet, wave-like gradient sequential inflation/deflation cycle provides patient with comfortable massage-like therapy.
- Easy-to-read LCD screen that allows the user to see real-time pressure in each individual chamber.
- Convenient digital controls, built-in timer feature for easy therapy time control.
- Automatic shut-off and emergency stop feature for maximum user safety assurance.
- Vacuum mechanism evacuates air from the garments upon shut-off of the system allowing for greater ease in the applying and taking off the garments.
- Quick-release hose connectors allow for easy removal of the garments for cleaning as needed.

### Specifications

**Dimension:** 10" (W) x 9.5" (D) x 4.5" (H)

**Weight:** 5.5 lbs

**Inflation Cycle:** 60 seconds

**Rest Time:** 12 seconds

**Cycle Time:** 15 seconds / chamber

**End Therapy Vacuum Time:** 180 seconds

**Treatment Time:** 0.5 – 3 Hours

**Pressure Range:** 20 – 80 mmHg

**Sleeve Prefill:** To 10 mmHg

**Gradient:** 7%; example 3 mmHg at 50 mmHg Distal Pressure

**Electrical:** 120VAC, 60Hz, 45 VA maximum

**Fuse Rated:** 250VAC, 1.0 AMP, SLO-BLO

**Protection Against Electric Shock:** Class II

**Operation Mode:** Intermittent Gradient Sequential