

CircuFlow™ Garment for Intermittent Pneumatic Compression (IPC) Pumps

Devon Medical Products CircuFlow™ Garment pairs with our CircuFlow™ 5100 and 5200 Intermittent Pneumatic Compression (IPC) Pumps. Garments are available for legs and for arms.

The CircuFlow™ models continuously pump and releases air into and out of the garment's four air chambers to create pressure that improves peripheral circulation and extracellular cleansing. The garment's design promotes high patient comfort and satisfaction.

Key Features

- Easy to apply and take off with easy-pull handles and mechanism that prevents full unzipping of the garment
- Lightweight, durable material for long-lasting use
- Easy to clean with instructions included to ensure proper care
- Ridging proof garment segment design that prevents tissue from being compressed in certain gap areas
- Reverse vacuum function that ensures easy removal of the garment after therapy



Model D-300S	Small leg garment (4 chambers)
Model D-300SW	Small wide leg garment (4 chambers)
Model D-300SXW	Small extra wide leg garment (4 chambers)
Model D-300M	Medium leg garment (4 chambers)
Model D-300MW	Medium wide leg garment (4 chambers)
Model D-300MXW	Medium extra wide leg garment (4 chambers)
Model D-300L	Large leg garment (4 chambers)
Model D-300LW	Large wide leg garment (4 chambers)
Model D-301R	Regular half leg garment (4 chambers)
Model D-301XL	Extra large half leg garment (4 chambers)
Model D-302M	Medium arm garment (4 chambers)
Model D-302L	Large arm garment (4 chambers)